

Opaguard™

High Level Instrument Disinfection

About Opaguard

Opaguard Solution is a convenient, ready-to-use, high-level disinfectant that provides a broad spectrum of kill. It is trusted by hospitals all the world over for providing rapid, cost-effective high-level disinfection for a wide range of endoscopes and other healthcare instruments.

Features & Benefits

- Glutaraldehyde-free (0.55% ortho-phthalaldehyde) high-level disinfecting solution.
- Rapid 5-minute soak time at 25°C in an automated endoscope reprocessor.
- Twelve-minute soak time 20°C for manual reprocessing.
- No activating or mixing needed.
- Minimizes exposure to fumes, odor or irritation
- Long-lasting efficacy allows reprocessing of more devices per gallon than with Glutaraldehyde.
- Effective against HIV, HBV, MTB & Spores.



Direction for use:

- Immerse device completely, filling all lumens and eliminating air pockets in Opaguard solution for a minimum of 5 minutes at 20°C or higher remove device from the solution and rinse thoroughly with sterile water or with potable tap water. Does not require activation before use.

Precautions:

- Use in well ventilated areas.
- Avoid direct skin contact.
- Pre-clean instruments with instruments cleanser after rinse under running water for optimum result.

Composition

Ortho-Phthalaldehyde 0.55% w/w

Technical Specification

pH	: 7.5 (Neutral)
Stability Temperature	: 15°C-30°C
Shelf life	: 2 years
Re-use life	: 14 days
Pack Size	: 1 Litre, 5 Litres

Genex Hygiene

B-104, First Floor Near, Munisurvart Compound, Rahnal Village, Bhiwandi, Thane
genexhygiene@gmail.com | www.genexhygiene.com