

# Instrugen™-2.45%

Sterilizing & Disinfecting Solution

## About Instrugen-2.45%

Instrugen-2.45% offers excellent material compatibility and can be used to disinfect a wide range of medical instruments, made of aluminum, brass, copper, stainless steel, plastics and elastomers.

Activated Instrugen-2.45% Solution provides a wide spectrum efficacy against bacteria, mycobacteria, viruses and fungi. The solution can also achieve some sporicidal activity with longer exposure time.



- **Rapid Action**
- **No Rusting**
- **No Staining**

- Feature & Benefits**
- Low purification index values represent higher purity
  - Ensures high purity, efficacy & stability
  - In accordance with CDC recommendation
  - Excellent material compatibility
  - Surfactant free
  - Enhanced use life

## Activity

### Broad spectrum:

- Bactericidal, Fungicidal, Virucidal & Sporicidal

### Contact Time

- High Level Disinfection** - 20 Minutes at 20°C
- Sporicidal Activity** - 10 Hours at 25°C
- Use Life** - Up to 28 days (Activated)
- Shelf Life** - 24 months
- Storage** - 15-30°C

## Composition:

- Glutaraldehyde : 2.45% w/v
- Inert Ingredients : 97.55%

## Applications

- GI Endoscopes
- Anaesthesia Equipments
- Diaphragm Fitting Rings
- Dental Burrs
- Hand Pieces

## Genex Hygiene

B-104, First Floor Near, Munisurvat Compound, Rahnal Village, Bhiwandi, Thane  
genexhygiene@gmail.com | www.genexhygiene.com